

For Research Use Only

# TTYH1 Recombinant antibody

Catalog Number: 80263-1-RR



## Basic Information

<b>Catalog Number:</b> 80263-1-RR	<b>GenBank Accession Number:</b> BC011347	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ul , Concentration: 1000 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 57348	<b>CloneNo.:</b> 8G21
<b>Source:</b> Rabbit	<b>Full Name:</b> tweety homolog 1 (Drosophila)	<b>Recommended Dilutions:</b> WB 1:5000-1:50000
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 49-50 kDa	
<b>Immunogen Catalog Number:</b> AG25290		

## Applications

**Tested Applications:**  
WB, ELISA

**Species Specificity:**  
Human, rabbit, pig

### Positive Controls:

WB : pig cerebellum tissue, human testis tissue, rabbit cerebellum tissue, rabbit brain tissue

## Background Information

### Storage


**Storage:**  
Store at -20°C.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 80263-1-RR (TTYH1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.