

For Research Use Only

# TTYH1 Recombinant antibody

Catalog Number: 80263-1-RR



## Basic Information

Catalog Number:

80263-1-RR

Size:

100ul , Concentration: 1000 ug/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25290

GenBank Accession Number:

BC011347

GeneID (NCBI):

57348

UNIPROT ID:

Q9H313

Full Name:

tweety homolog 1 (Drosophila)

Calculated MW:

49-50 kDa

Purification Method:

Protein A purification

CloneNo.:

8G21

Recommended Dilutions:

WB 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, rabbit, pig

Positive Controls:

WB : pig cerebellum tissue, human testis tissue, rabbit cerebellum tissue, rabbit brain tissue

## Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 80263-1-RR (TTYH1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.