

For Research Use Only

PIK3CG Recombinant monoclonal antibody

Catalog Number: 80720-2-RR



Basic Information

Catalog Number:

80720-2-RR

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002649

GeneID (NCBI):

5294

UNIPROT ID:

P48736

Full Name:

phosphoinositide-3-kinase, catalytic, gamma polypeptide

Calculated MW:

126 kDa

Observed MW:

100-130 kDa

Purification Method:

Protein A purification

CloneNo.:

251747C3

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : K-562 cells,

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 80720-2-RR (PIK3CG antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.