

For Research Use Only



DYKDDDDK tag Recombinant antibody, PBS Only (Binds to FLAG® tag epitope)

Catalog Number: **80745-1-PBS**

Basic Information

Catalog Number:

80745-1-PBS

Size:

100ug , Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2329

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

5E11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

recombinant protein

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



TDP-43-flag Transfected HEK-293T cells were subjected to SDS PAGE followed by western blot with 80745-1-RR (DYKDDDDK tag antibody) at dilution of 1:50000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 80745-1-PBS in a different storage buffer formulation.