

For Research Use Only

XBP1 Recombinant antibody

Catalog Number: 80747-1-RR

1 Publications



Basic Information

Catalog Number:

80747-1-RR

Size:

100ul

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21454

GenBank Accession Number:

BC000938

GeneID (NCBI):

7494

UNIPROT ID:

P17861

Full Name:

X-box binding protein 1

Calculated MW:

261 aa, 29 kDa

Purification Method:

Protein A purification

CloneNo.:

117

Recommended Dilutions:

WB 1:2000-1:18000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human, mouse

Positive Controls:

WB : HEK-293 cells, Jurkat cells

Notable Publications

Author	Pubmed ID	Journal	Application
Zhou Du	35986482	J Pineal Res	WB

Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 80747-1-RR (XBP1 antibody) at dilution of 1:9000 incubated at room temperature for 1.5 hours.