

For Research Use Only

# FGFR4 Recombinant antibody

Catalog Number: 81070-1-RR



## Basic Information

<b>Catalog Number:</b> 81070-1-RR	<b>GenBank Accession Number:</b> BC011847	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ul , Concentration: 2292 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 2264	<b>CloneNo.:</b> 2M22
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P22455	<b>Recommended Dilutions:</b> WB 1:4000-1:40000
<b>Isotype:</b> IgG	<b>Full Name:</b> fibroblast growth factor receptor 4	
<b>Immunogen Catalog Number:</b> AG31013	<b>Calculated MW:</b> 88 kDa	
	<b>Observed MW:</b> 100-110 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b>
<b>Species Specificity:</b> Human	<b>WB :</b> MOLT-4 cells, Raji cells, Jurkat cells, HeLa cells, HepG2 cells, HT-29 cells, A549 cells

## Storage


**Storage:**  
Store at -20°C.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 81070-1-RR (FGFR4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.