

For Research Use Only

FASN Recombinant monoclonal antibody, PBS Only

Catalog Number: 81079-2-PBS



Basic Information

Catalog Number:

81079-2-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0975

GenBank Accession Number:

BC007909

GeneID (NCBI):

2194

UNIPROT ID:

P49327

Full Name:

fatty acid synthase

Calculated MW:

272 kDa

Observed MW:

250-272 kDa

Purification Method:

Protein A purification

CloneNo.:

260308B2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

FASN gene codes for an enzyme essential for de novo fatty acid synthesis and cellular substrate energy metabolism. Active FASN is a homodimer in which each peptide subunit has a molecular weight of 260 kDa. FASN is overexpressed in various types of cancer including glioblastomas and is a potential therapeutic target. Recently FASN has been reported to contribute to the neurogenesis since FASN mutation caused intellectual disability in mice.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 81079-2-RR (FASN antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 81079-2-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 81079-2-RR (FASN antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 81079-2-PBS in a different storage buffer formulation.



rat liver tissue were subjected to SDS PAGE followed by western blot with 81079-2-RR (FASN antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 81079-2-PBS in a different storage buffer formulation.