

For Research Use Only

RAB1A Recombinant antibody

Catalog Number: 81837-1-RR



Basic Information

Catalog Number: 81837-1-RR	GenBank Accession Number: BC000905	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 5861	CloneNo.: 2G24
Source: Rabbit	Full Name: RAB1A, member RAS oncogene family	Recommended Dilutions: WB 1:5000-1:50000
Isotype: IgG	Calculated MW: 205 aa, 23 kDa	
Immunogen Catalog Number: AG2264	Observed MW: 23-25 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human, mouse, rat

Positive Controls:

WB : HeLa cells, HepG2 cells, mouse brain tissue, rat brain tissue

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 81837-1-RR (RAB1A antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.