

For Research Use Only

CHAT Recombinant antibody, PBS Only (Capture)



Catalog Number: 81994-2-PBS

Basic Information

Catalog Number: 81994-2-PBS	GenBank Accession Number: NM_020549	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 1103	CloneNo.: 240604E1
Source: Rabbit	UNIPROT ID: P28329	
Isotype: IgG	Full Name: choline acetyltransferase	
	Calculated MW: 83 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

81994-2-PBS targets CHAT as part of a matched antibody pair.

MP00576-2: 81994-2-PBS capture and 81994-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00576-2, CHAT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81994-2-PBS. Detection antibody: 81994-4-PBS. Standard: SY00995. Range: 19.5-2500 pg/mL.