

For Research Use Only

Desmin Recombinant monoclonal antibody, PBS Only

Catalog Number: 82748-2-PBS



Basic Information

Catalog Number: 82748-2-PBS	GenBank Accession Number: BC032116	Purification Method: Protein A purification
Size: 100ug, Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 1674	CloneNo.: 5K5
Source: Rabbit	UNIPROT ID: P17661	
Isotype: IgG	Full Name: desmin	
Immunogen Catalog Number: AG9742	Calculated MW: 470 aa, 54 kDa	

Applications

Tested Applications:
IHC, Indirect ELISA

Species Specificity:
human

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 82748-2-RR (Desmin antibody) at dilution of 1:3000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82748-2-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 82748-2-RR (Desmin antibody) at dilution of 1:3000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82748-2-PBS in a different storage buffer formulation.