

For Research Use Only

MUC2 Recombinant monoclonal antibody, PBS Only

Catalog Number: 82792-2-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 82792-2-PBS | GenBank Accession Number: M94132 | Purification Method: Protein A purification |
| Size: 100ug , Concentration: 1 mg/ml by Nanodrop; | GeneID (NCBI): 4583 | CloneNo.: 250669F4 |
| Source: Rabbit | Full Name: mucin 2, oligomeric mucus/gel-forming | |
| Isotype: IgG | Calculated MW: 540 kDa | |
| Immunogen Catalog Number: AG25800 | | |

Applications

Tested Applications:
IHC, IF-P, Indirect ELISA

Species Specificity:
human, mouse

Background Information

This gene encodes a member of the mucin protein family. Mucins are high molecular weight glycoproteins produced by many epithelial tissues. The protein encoded by this gene is secreted and forms an insoluble mucous barrier that protects the gut lumen. The protein polymerizes into a gel of which 80% is composed of oligosaccharide side chains by weight. Downregulation of this gene has been observed in patients with Crohn disease and ulcerative colitis.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 82792-2-RR (MUC2 antibody) at dilution of 1:2000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82792-2-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 82792-2-RR (MUC2 antibody) at dilution of 1:2000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82792-2-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using 82792-2-RR (MUC2 antibody) at dilution of 1:2000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82792-2-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using 82792-2-RR (MUC2 antibody) at dilution of 1:2000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82792-2-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 82792-2-RR (MUC2 antibody) at dilution of 1:2000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82792-2-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded human colon tissue using MUC2 antibody (82792-2-RR, Clone: 250669F4) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82792-2-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse colon tissue using MUC2 antibody (82792-2-RR, Clone: 250669F4) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82792-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82792-2-RR against Human MUC2 were performed. The affinity constant is below 1 pM.