

For Research Use Only

# TNFRSF17 Recombinant antibody

Catalog Number: 82855-1-RR



## Basic Information

<b>Catalog Number:</b> 82855-1-RR	<b>GenBank Accession Number:</b> BC058291	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ul , Concentration: 300 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 608	<b>CloneNo.:</b> 5H7
<b>Source:</b> Rabbit	<b>Full Name:</b> tumor necrosis factor receptor superfamily, member 17	<b>Recommended Dilutions:</b> WB 1:1000-1:4000
<b>Isotype:</b> IgG	<b>Calculated MW:</b> 184 aa, 20 kDa	
	<b>Observed MW:</b> 20 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : Ramos cells, IM-9 cells
<b>Species Specificity:</b> Human	

## Background Information

### Storage


**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

**\*\*\* 20ul sizes contain 0.1%BSA**

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)      E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82855-1-RR (TNFRSF17 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.