

For Research Use Only

Mouse IL-12/IL-23 p40 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82865-9-PBS



Basic Information

Catalog Number:

82865-9-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001303244

GeneID (NCBI):

16160

UNIPROT ID:

P43432

Full Name:

interleukin 12b

Calculated MW:

38kd

Purification Method:

Protein A purification

CloneNo.:

240210H6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Product Information

82865-9-PBS targets IL-12/IL-23 p40 as part of a matched antibody pair:

MP00413-2: 82865-9-PBS capture and 82865-7-PBS detection (validated in Cytometric bead array)

MP00413-3: 82865-9-PBS capture and 82865-6-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00413-2, MOUSE IL-12/IL-23 p40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82865-9-PBS. Detection antibody: 82865-7-PBS. Standard: Eg0299. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00413-3, MOUSE IL-12/IL-23 p40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82865-9-PBS. Detection antibody: 82865-6-PBS. Standard: Eg0299. Range: 5-80 ng/mL