

For Research Use Only

# SFTPb Recombinant antibody

Catalog Number: 82866-16-RR



## Basic Information

Catalog Number:

82866-16-RR

Size:

100ul, Concentration: 250 ug/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23605

GenBank Accession Number:

BC032785

GeneID (NCBI):

6439

UNIPROT ID:

P07988

Full Name:

surfactant protein B

Calculated MW:

381 aa, 42 kDa

Observed MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

230240D4

Recommended Dilutions:

WB 1:1000-1:5000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB : mouse lung tissue,

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse lung tissue were subjected to SDS PAGE followed by western blot with 82866-16-RR (SFTPB antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.