

For Research Use Only

SFTPB Recombinant antibody

Catalog Number: 82866-5-RR



Basic Information

Catalog Number: 82866-5-RR	GenBank Accession Number: BC032785	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 6439	CloneNo.: 230240E5
Source: Rabbit	UNIPROT ID: P07988	
Isotype: IgG	Full Name: surfactant protein B	
Immunogen Catalog Number: AG23605	Calculated MW: 381 aa, 42 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

