

For Research Use Only

SFTPB Recombinant antibody

Catalog Number: 82866-5-RR



Basic Information

Catalog Number: 82866-5-RR	GenBank Accession Number: BC032785	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 6439	CloneNo.: 230240E5
Source: Rabbit	Full Name: surfactant protein B	
Isotype: IgG	Calculated MW: 381 aa, 42 kDa	
Immunogen Catalog Number: AG23605		

Applications

Tested Applications:
FC, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ A549 cells were intracellularly stained with 0.13 ug Anti-Human SFTPB (82866-5-RR, Clone:230240E5) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) at dilution 1:166 (red), or 0.13 ug Rabbit IgG control Rabbit PolyAb (30000-0-AP, Clone:) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).