

For Research Use Only

CHD5 Recombinant antibody

Catalog Number: 82896-1-RR



Basic Information

Catalog Number: 82896-1-RR	GenBank Accession Number: BC148803	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 26038	CloneNo.: 230093A6
Source: Rabbit	Full Name: chromodomain helicase DNA binding protein 5	
Isotype: IgG	Calculated MW: 1954 aa, 223 kDa	
Immunogen Catalog Number: AG18621		

Applications

Tested Applications:
FC, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.


Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ MCF-7 cells were intracellularly stained with 0.13 ug Anti-Human CHD5 (82896-1-RR, Clone:230093A6) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) at dilution 1:166 (red), or 0.13 ug Rabbit IgG control Rabbit PolyAb (30000-0-AP, Clone:) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).