

For Research Use Only

# CHD5 Recombinant antibody

Catalog Number: 82896-1-RR



## Basic Information

|  |  |   |
|--|--|---|
| <b>Catalog Number:</b><br>82896-1-RR                           | <b>GenBank Accession Number:</b><br>BC148803                     | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>100ul , Concentration: 1000 µg/ml by Nanodrop; | <b>GeneID (NCBI):</b><br>26038                                   | <b>CloneNo.:</b><br>230093A6                          |
| <b>Source:</b><br>Rabbit                                       | <b>UNIPROT ID:</b><br>Q8TD10                                     |   |
| <b>Isotype:</b><br>IgG   | <b>Full Name:</b><br>chromodomain helicase DNA binding protein 5 |   |
| <b>Immunogen Catalog Number:</b><br>AG18621                    | <b>Calculated MW:</b><br>1954 aa, 223 kDa                        |   |

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
Human

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

