

For Research Use Only

ALPL Recombinant antibody

Catalog Number: 82915-1-RR



Basic Information

Catalog Number: 82915-1-RR	GenBank Accession Number: BC021289	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 249	CloneNo.: 230212B10
Source: Rabbit	UNIPROT ID: P05186	
Isotype: IgG	Full Name: alkaline phosphatase, liver/bone/kidney	
Immunogen Catalog Number: AG1664	Calculated MW: 57 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

