

For Research Use Only

MYOF Recombinant antibody, PBS Only



Catalog Number: 82917-2-PBS

Basic Information

Catalog Number: 82917-2-PBS	GenBank Accession Number: NM_133337	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 26509	CloneNo.: 230213E10
Source: Rabbit	UNIPROT ID: Q9NZM1	
Isotype: IgG	Full Name: myoferlin	
	Calculated MW: 235 kDa	
	Observed MW: 230-240 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human

Background Information

MYOF, also named as FER1L3 and KIAA1207, belongs to the ferlin family. MYOF is a calcium/phospholipid-binding protein that plays a role in the plasmalemma repair mechanism of endothelial cells that permits rapid resealing of membranes disrupted by mechanical stress. It is involved in endocytic recycling. MYOF is implicated in VEGF signal transduction by regulating the levels of the receptor KDR.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82917-2-RR (Myoferlin antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82917-2-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using MYOF antibody (82917-2-RR, Clone: 230213E10) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 82917-2-PBS in a different storage buffer formulation.



1×10^6 HeLa cells were intracellularly stained with 0.25 ug MYOF Recombinant antibody (82917-2-RR, Clone:230213E10) and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2) (red), or 0.25 ug Isotype Control (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 82917-2-PBS in a different storage buffer formulation.