

For Research Use Only

MYOF Recombinant antibody

Catalog Number: 82917-7-RR



Basic Information

Catalog Number:

82917-7-RR

GenBank Accession Number:

NM_133337

Purification Method:

Protein A purification

Size:

100ul , Concentration: 1000 ug/ml by

GeneID (NCBI):

26509

CloneNo.:

230213B6

Nanodrop;

UNIPROT ID:

Q9NZM1

Recommended Dilutions:

IF/ICC 1:125-1:500

Source:

Rabbit

Full Name:

myoferlin

Isotype:

IgG

Calculated MW:

235 kDa

Applications

Tested Applications:

IF/ICC, ELISA

Positive Controls:

IF/ICC : HeLa cells,

Species Specificity:

Human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using MYOF antibody (82917-7-RR, Clone: 230213B6) at dilution of 1:250 and CoraLite®488-Conjugated Goat Anti-Mouse IgG(H+L) (SA00013-1).