

For Research Use Only

# MAPKAPK2 Recombinant antibody

Catalog Number: 82929-2-RR



## Basic Information

<b>Catalog Number:</b> 82929-2-RR	<b>GenBank Accession Number:</b> BC036060	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ul , Concentration: 700 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 9261	<b>CloneNo.:</b> 230197D6
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P49137	<b>Recommended Dilutions:</b> WB 1:2000-1:14000
<b>Isotype:</b> IgG	<b>Full Name:</b> mitogen-activated protein kinase-activated protein kinase 2	
<b>Immunogen Catalog Number:</b> AG34196	<b>Calculated MW:</b> 400 aa, 46 kDa	
	<b>Observed MW:</b> 42-46 kDa	

## Applications

<b>Tested Applications:</b> WB, ELISA	<b>Positive Controls:</b> WB : A549 cells, HeLa cells, HEK-293T cells
<b>Species Specificity:</b> Human	

## Storage


<b>Storage:</b> Store at -20°C. Stable for one year after shipment.
<b>Storage Buffer:</b> PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82929-2-RR (MAPKAPK2 antibody) at dilution of 1:7000 incubated at room temperature for 1.5 hours.