

For Research Use Only

RBP4 Recombinant antibody

Catalog Number: 82934-1-RR



Basic Information

Catalog Number: 82934-1-RR	GenBank Accession Number: BC020633	Purification Method: Protein A purification
Size: 100ul , Concentration: 400 µg/ml by Nanodrop;	GeneID (NCBI): 5950	CloneNo.: 6013
Source: Rabbit	Full Name: retinol binding protein 4, plasma	Recommended Dilutions: WB 1:1000-1:4000
Isotype: IgG	Calculated MW: 201 aa, 23 kDa	
Immunogen Catalog Number: AG2448	Observed MW: 23 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : HepG2 cells,
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82934-1-RR (RBP4 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.