

For Research Use Only

TAT Recombinant antibody, PBS Only

Catalog Number: 82939-1-PBS



Basic Information

Catalog Number:

82939-1-PBS

Size:

100ug , Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34043

GenBank Accession Number:

BC130534

GeneID (NCBI):

6898

UNIPROT ID:

P17735

Full Name:

tyrosine aminotransferase

Observed MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

230314H2

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

Human, mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

