

For Research Use Only

TAT Recombinant antibody

Catalog Number: 82939-1-RR



Basic Information

Catalog Number: 82939-1-RR	GenBank Accession Number: BC130534	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 6898	CloneNo.: 230314H2
Source: Rabbit	Full Name: tyrosine aminotransferase	Recommended Dilutions: WB 1:2000-1:10000
Isotype: IgG	Observed MW: 60 kDa	
Immunogen Catalog Number: AG34043		

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Raji cells,
Species Specificity: Human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82939-1-RR (TAT antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.