

For Research Use Only

SCN3B Recombinant antibody, PBS Only



Catalog Number: 82959-7-PBS

Basic Information

Catalog Number:

82959-7-PBS

Size:

100ug , Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23201

GenBank Accession Number:

BC126265

GeneID (NCBI):

55800

UNIPROT ID:

Q9NY72

Full Name:

sodium channel, voltage-gated, type III, beta

Purification Method:

Protein A purification

CloneNo.:

230316B2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82959-7-RR (SCN3B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82959-7-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82959-7-RR against Human SCN3B were performed. The affinity constant is 0.98 nM.