

For Research Use Only

ACAP1 Recombinant antibody, PBS Only

Catalog Number: 82967-1-PBS

Featured Product



Basic Information

Catalog Number:

82967-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34309

GenBank Accession Number:

BC018543

GeneID (NCBI):

9744

UNIPROT ID:

Q15027

Full Name:

ArfGAP with coiled-coil, ankyrin repeat and PH domains 1

Calculated MW:

740 aa, 82 kDa

Observed MW:

82-85 kDa

Purification Method:

Protein A purification

CloneNo.:

230370E1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

ACAP1, also named as Centaurin-beta-1, is a 740 amino acid protein, which is expressed Highly in lung and spleen. ACAP1 as a GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6) is required for clathrin-dependent export of proteins from recycling endosomes to trans-Golgi network and cell surface and required for regulated export of ITGB1 from recycling endosomes to the cell surface and ITGB1-dependent cell migration.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




WB result of ACAP1 antibody (82967-1-RR; 1:3000; incubated at room temperature for 1.5 hours) with sh-Control and sh-ACAP1 transfected Jurkat cells. This data was developed using the same antibody clone with 82967-1-PBS in a different storage buffer formulation.



Jurkat cells were subjected to SDS PAGE followed by western blot with 82967-1-RR (ACAP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82967-1-PBS in a different storage buffer formulation.



K-562 cells were subjected to SDS PAGE followed by western blot with 82967-1-RR (ACAP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82967-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82967-1-RR against Human ACAP1 were performed. The affinity constant is 0.42 nM.