

ARID1B Recombinant antibody, PBS Only (Detector)

Catalog Number: 82979-5-PBS

Basic Information

Catalog Number: 82979-5-PBS	GenBank Accession Number: NM_017519	Purification Method: Protein A purification
Size: 100ug, Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 57492	CloneNo.: 230222D7
Source: Rabbit	UNIPROT ID: Q8NFD5	
Isotype: IgG	Full Name: AT rich interactive domain 1B (SWI1-like)	
Immunogen Catalog Number: AG30643	Calculated MW: 236 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

82979-5-PBS targets ARID1B as part of a matched antibody pair:

MP00103-1: 82979-4-PBS capture and 82979-5-PBS detection (validated in Cytometric bead array)

MP00103-2: 82979-2-PBS capture and 82979-5-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00103-1, ARID1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82979-4-PBS. Detection antibody: 82979-5-PBS. Standard: Ag30643. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00103-2, ARID1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82979-2-PBS. Detection antibody: 82979-5-PBS. Standard: Ag30643. Range: 0.625-20 ng/mL