

For Research Use Only

ADAMDEC1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 82985-5-PBS



Basic Information

Catalog Number: 82985-5-PBS	GenBank Accession Number: BC074877	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 27299	CloneNo.: 230367D12
Source: Rabbit	UNIPROT ID: O15204	
Isotype: IgG	Full Name: ADAM-like, decysin 1	
Immunogen Catalog Number: AG12265	Calculated MW: 470 aa, 53 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

82985-5-PBS targets ADAMDEC1 as part of a matched antibody pair:

MP00007-1: 82985-5-PBS capture and 82985-4-PBS detection (validated in Cytometric bead array)

MP00007-2: 82985-5-PBS capture and 82985-3-PBS detection (validated in Cytometric bead array)

MP00007-3: 82985-2-PBS capture and 82985-5-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00007-2, ADAMDEC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82985-5-PBS. Detection antibody: 82985-3-PBS. Standard:Ag12265. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00007-1, ADAMDEC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82985-5-PBS. Detection antibody: 82985-4-PBS. Standard:Ag12265. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00007-3, ADAMDEC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82985-2-PBS. Detection antibody: 82985-5-PBS. Standard:Ag12265. Range: 0.625-20 ng/mL