

For Research Use Only

Mouse Hif1a Recombinant antibody, PBS Only

Catalog Number: 82989-4-PBS



Basic Information

Catalog Number:

82989-4-PBS

Size:

100ug, Concentration: 1mg/ml by
Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33620

GenBank Accession Number:

NM_001313919

GeneID (NCBI):

15251

UNIPROT ID:

Q61221

Full Name:

hypoxia inducible factor 1, alpha
subunit

Calculated MW:

93 kDa

Purification Method:

Protein A purification

CloneNo.:

230101E8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Mouse

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Cobalt Chloride treated NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 82989-4-RR (Hif1a antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82989-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82989-4-RR against Mouse Hif1a were performed. The affinity constant is below 1 pM.