

For Research Use Only

POLR2K Recombinant antibody, PBS Only (Detector)

Catalog Number: 82995-1-PBS



Basic Information

Catalog Number: 82995-1-PBS	GenBank Accession Number: BC018157	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 5440	CloneNo.: 230147E8
Source: Rabbit	UNIPROT ID: P53803	
Isotype: IgG	Full Name: polymerase (RNA) II (DNA directed) polypeptide K, 7.0kDa	
Immunogen Catalog Number: AG33765		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

82995-1-PBS targets POLR2K as part of a matched antibody pair:

MP01314-1: 82995-3-PBS capture and 82995-1-PBS detection (validated in Cytometric bead array)

MP01314-2: 82995-2-PBS capture and 82995-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01314-1, POLR2K Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82995-3-PBS. Detection antibody: 82995-1-PBS. Standard: Ag33765. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP01314-2, POLR2K Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82995-2-PBS. Detection antibody: 82995-1-PBS. Standard: Ag33765. Range: 0.313-20 ng/mL