

For Research Use Only

# UNC13A Recombinant antibody, PBS Only



Catalog Number: 83000-6-PBS

## Basic Information

<b>Catalog Number:</b> 83000-6-PBS	<b>GenBank Accession Number:</b> NM_001080421	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug, Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 23025	<b>CloneNo.:</b> 230469G2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9UPW8	
<b>Isotype:</b> IgG	<b>Full Name:</b> unc-13 homolog A (C. elegans)	
	<b>Calculated MW:</b> 193 kDa	

## Applications

**Tested Applications:**  
IF/ICC, ELISA

**Species Specificity:**  
human

## Storage


**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed SH-SY5Y cells using UNC13A/Munc13-1 antibody (83000-6-RR, Clone: 230469G2 ) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 83000-6-PBS in a different storage buffer formulation.