

For Research Use Only

GCSH Recombinant antibody, PBS Only



Catalog Number: 83012-5-PBS

Basic Information

Catalog Number: 83012-5-PBS	GenBank Accession Number: BC000790	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 2653	CloneNo.: 230505A3
Source: Rabbit	UNIPROT ID: P23434	
Isotype: IgG	Full Name: glycine cleavage system protein H (aminomethyl carrier)	
Immunogen Catalog Number: AG10174	Calculated MW: 19 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 83012-5-RR (GCSH antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83012-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83012-5-RR against Human GCSH were performed. The affinity constant is below 1 pM.