

For Research Use Only

# Cadherin-19 Recombinant antibody

Catalog Number: 83019-1-RR



## Basic Information

<b>Catalog Number:</b> 83019-1-RR	<b>GenBank Accession Number:</b> BC015877	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ul , Concentration: 1000 µg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 28513	<b>CloneNo.:</b> 230271D2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9H159	
<b>Isotype:</b> IgG	<b>Full Name:</b> cadherin 19, type 2	
<b>Immunogen Catalog Number:</b> AG2000	<b>Calculated MW:</b> 55 kDa; 87 kDa	

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
Human

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

