

For Research Use Only

BMPR1B Recombinant antibody, PBS Only (Detector)

Catalog Number: 83029-1-PBS



Basic Information

Catalog Number: 83029-1-PBS	GenBank Accession Number: BC047773	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 658	CloneNo.: 230244B11
Source: Rabbit	UNIPROT ID: O00238	Affinity: $K_D = 1.49 \times 10^{-10}M$
Isotype: IgG	Full Name: bone morphogenetic protein receptor, type IB	$K_{off} = 1.22 \times 10^{-4}M$ $K_{on} = 8.14 \times 10^4M$
Immunogen Catalog Number: AG24902	Calculated MW: 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

83029-1-PBS targets BMPR1B as part of a matched antibody pair:

MP00108-2: 83029-3-PBS capture and 83029-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00108-2, BMPRII Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83029-3-PBS. Detection antibody: 83029-1-PBS. Standard: Ag24902. Range: 0.313-10 ng/mL.