

For Research Use Only

SLC29A3 Recombinant antibody

Catalog Number: 83036-3-RR



Basic Information

Catalog Number:

83036-3-RR

Size:

100ul, Concentration: 1000 µg/ml by

Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18563

GenBank Accession Number:

BC120996

GeneID (NCBI):

55315

UNIPROT ID:

Q9BZD2

Full Name:

solute carrier family 29 (nucleoside transporters), member 3

Calculated MW:

475 aa, 52 kDa

Observed MW:

51 kDa, 60 kDa

Purification Method:

Protein A purification

CloneNo.:

240123A8

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, A549 cells, LNCaP cells

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83036-3-RR (SLC29A3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 83036-3-RR against Human SLC29A3 were performed. The affinity constant is 0.824 nM.