

For Research Use Only

KIAA0776 Recombinant antibody

Catalog Number: 83047-6-RR



Basic Information

Catalog Number:

83047-6-RR

GenBank Accession Number:

BC036379

Purification Method:

Protein A purification

Size:

100ul , Concentration: 1000 µg/ml by Nanodrop;

GeneID (NCBI):

23376

CloneNo.:

230363A8

Source:

Rabbit

UNIPROT ID:

O94874

Isotype:

IgG

Full Name:

KIAA0776

Immunogen Catalog Number:

AG23362

Applications

Tested Applications:

ELISA

Species Specificity:

Human

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

