

For Research Use Only

# Integrin beta 3 Recombinant antibody

Catalog Number: 83053-5-RR

1 Publications



## Basic Information

Catalog Number:

83053-5-RR

Size:

100ul

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC127667

GeneID (NCBI):

3690

UNIPROT ID:

P05106

Full Name:

integrin, beta 3 (platelet glycoprotein IIIa, antigen CD61)

Calculated MW:

87 kDa

Purification Method:

Protein A purification

CloneNo.:

230292H2

## Applications

Tested Applications:

ELISA

Cited Applications:

WB

Species Specificity:

Human

Cited Species:

human

## Notable Publications

Author	Pubmed ID	Journal	Application
Jingyuan Zhang	39754146	Chin Med	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

\*\*\* 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

