

For Research Use Only

FEZF1 Recombinant monoclonal antibody, PBS Only

Catalog Number: 83061-1-PBS



Basic Information

Catalog Number:

83061-1-PBS

Size:

100ug, Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001024613

GeneID (NCBI):

389549

UNIPROT ID:

A0PJY2

Full Name:

FEZ family zinc finger 1

Calculated MW:

52 kDa

Purification Method:

Protein A purification

CloneNo.:

230341E1

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

FEZF1 is a highly conserved transcription factor belonging to the large family of C2H2 zinc finger proteins. FEZF1 is expressed early during mouse development and is important for brain development and cell identity.

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data