

For Research Use Only

FEZF1 Recombinant monoclonal antibody

Catalog Number: 83061-1-RR



Basic Information

| | | |
|----------------------------------------------------------------|--------------------------------------------------|-------------------------------------------------------|
| Catalog Number: 83061-1-RR | GenBank Accession Number: NM_001024613 | Purification Method: Protein A purification |
| Size: 100ul , Concentration: 1000 µg/ml by Nanodrop; | GeneID (NCBI): 389549 | CloneNo.: 230341E1 |
| Source: Rabbit | UNIPROT ID: A0PJY2 | |
| Isotype: IgG | Full Name: FEZ family zinc finger 1 | |
| | Calculated MW: 52 kDa | |

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

FEZF1 is a highly conserved transcription factor belonging to the large family of C2H2 zinc finger proteins. FEZF1 is expressed early during mouse development and is important for brain development and cell identity.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

