

For Research Use Only

STT3A Recombinant antibody, PBS Only



Catalog Number: 83077-4-PBS

Basic Information

Catalog Number:

83077-4-PBS

Size:

100ug, Concentration: 1mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31994

GenBank Accession Number:

BC020965

GeneID (NCBI):

3703

UNIPROT ID:

P46977

Full Name:

STT3, subunit of the oligosaccharyltransferase complex, homolog A (*S. cerevisiae*)

Calculated MW:

705 aa, 81 kDa

Purification Method:

Protein A purification

CloneNo.:

230426B10

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 83077-4-RR against Human STT3A were performed. The affinity constant is 9.73 nM.