

For Research Use Only

CELSR2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83115-3-PBS



Basic Information

Catalog Number: 83115-3-PBS	GenBank Accession Number: NM_001408	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 1952	CloneNo.: 230319D10
Source: Rabbit	UNIPROT ID: Q9HCU4	
Isotype: IgG	Full Name: cadherin, EGF LAG seven-pass G-type receptor 2 (flamingo homolog, Drosophila)	
	Calculated MW: 318 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Product Information

83115-3-PBS targets CELSR2 as part of a matched antibody pair:

MP00165-2: 83115-2-PBS capture and 83115-3-PBS detection (validated in Cytometric bead array, Sandwich ELISA)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00165-1, Human CELSR2
Recombinant Matched Antibody Pair, PBS Only
measured by Cytometric bead array. Capture
antibody: 83115-4-PBS. Detection antibody: 83115-
3-PBS. Standard:null. Range: null.