

For Research Use Only

# CELSR2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83115-4-PBS



## Basic Information

**Catalog Number:**

83115-4-PBS

**Size:**

100ug, Concentration: 1 mg/ml by Nanodrop;

**Source:**

Rabbit

**Isotype:**

IgG

**GenBank Accession Number:**

NM\_001408

**GeneID (NCBI):**

1952

**UNIPROT ID:**

Q9HCU4

**Full Name:**

cadherin, EGF LAG seven-pass G-type receptor 2 (flamingo homolog, Drosophila)

**Calculated MW:**

318 kDa

**Purification Method:**

Protein A purification

**CloneNo.:**

230319F4

## Applications

**Tested Applications:**

Cytometric bead array, Indirect ELISA

**Species Specificity:**

human

## Storage

**Storage:**

Store at -80°C.

**Storage Buffer:**

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00165-1, Human CELSR2  
Recombinant Matched Antibody Pair, PBS Only  
measured by Cytometric bead array. Capture  
antibody: 83115-4-PBS. Detection antibody: 83115-  
3-PBS. Standard:null. Range: null.