

For Research Use Only

NLRP10 Recombinant antibody, PBS Only

Catalog Number: 83131-5-PBS



Basic Information

Catalog Number: 83131-5-PBS	GenBank Accession Number: BC104957	Purification Method: Protein A purification
Size: 100ug , Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 338322	CloneNo.: 230402D9
Source: Rabbit	UNIPROT ID: Q86W26	
Isotype: IgG	Full Name: NLR family, pyrin domain containing 10	
Immunogen Catalog Number: AG26402	Calculated MW: 655 aa, 75 kDa	

Applications

Tested Applications:
Indirect ELISA

Species Specificity:
human

Product Information

83131-5-PBS targets NLRP10 as part of a matched antibody pair:

MP00242-1: 83131-5-PBS capture and 83131-6-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

