

GBX2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 83145-4-PBS

Basic Information

Catalog Number: 83145-4-PBS	GenBank Accession Number: BC137448	Purification Method: Protein A purification
Size: 100ug, Concentration: 1mg/ml by Nanodrop;	GeneID (NCBI): 2637	CloneNo.: 230468G8
Source: Rabbit	UNIPROT ID: P52951	
Isotype: IgG	Full Name: gastrulation brain homeobox 2	
Immunogen Catalog Number: AG16288	Calculated MW: 348 aa, 37 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Product Information

83145-4-PBS targets GBX2 as part of a matched antibody pair:

MP00200-1: 83145-2-PBS capture and 83145-4-PBS detection (validated in Cytometric bead array)

MP00200-2: 83145-3-PBS capture and 83145-4-PBS detection (validated in Cytometric bead array)

MP00200-3: 83145-4-PBS capture and 83145-1-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.


Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve ofMP00200-1, GBX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83145-2-PBS. Detection antibody: 83145-4-PBS. Standard: Ag16288. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00200-2, GBX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83145-3-PBS. Detection antibody: 83145-4-PBS. Standard: Ag16288. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00200-3, GBX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83145-4-PBS. Detection antibody: 83145-1-PBS. Standard: Ag16288. Range: 0.625-20 ng/mL