

For Research Use Only

IL-17RB Recombinant antibody

Catalog Number: 83149-3-RR



Basic Information

Catalog Number: 83149-3-RR	GenBank Accession Number: NM_018725	Purification Method: Protein A purification
Size: 100ul , Concentration: 1000 µg/ml by Nanodrop;	GeneID (NCBI): 55540	CloneNo.: 230439B3
Source: Rabbit	UNIPROT ID: Q9NRM6	
Isotype: IgG	Full Name: interleukin 17 receptor B	
	Calculated MW: 56 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Background Information

IL17RB, also named as CRL4, EVI27 and IL17BR, is a receptor for the proinflammatory cytokines IL17B and IL17E. IL17RB may play a role in controlling the growth and/or differentiation of hematopoietic cells. IL17RB has two isoforms with MW 56 kDa (isoform 1) and 32 kDa (isoform 2). The antibody is specific to IL17RB.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

*** 20ul sizes contain 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data