

For Research Use Only

# STX1A Recombinant antibody, PBS Only

Catalog Number: 83159-6-PBS



## Basic Information

Catalog Number:

83159-6-PBS

Size:

100ug , Concentration: 1 mg/ml by Nanodrop;

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_004603

GeneID (NCBI):

6804

UNIPROT ID:

Q16623

Full Name:

syntaxin 1A (brain)

Calculated MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

230470F5

## Applications

Tested Applications:

WB, IF-P, ELISA

Species Specificity:

rat, mouse, human

## Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83159-6-RR (STX1A antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83159-6-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed mouse cerebellum tissue using STX1A antibody (83159-6-RR, Clone: 230470F5) at dilution of 1:500 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 83159-6-PBS in a different storage buffer formulation.