

For Research Use Only

COX6B2 Recombinant antibody, PBS Only



Catalog Number: 83164-3-PBS

Basic Information

Catalog Number: 83164-3-PBS	GenBank Accession Number: BC064548	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 125965	CloneNo.: 230497H6
Source: Rabbit	UNIPROT ID: Q6YFQ2	
Isotype: IgG	Full Name: cytochrome c oxidase subunit VIb polypeptide 2 (testis)	
Immunogen Catalog Number: AG1983	Calculated MW: 88 aa, 11 kDa	

Applications

Tested Applications:
ELISA, WB

Species Specificity:
mouse, human

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 83164-3-RR (COX6B2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83164-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83164-3-RR against Human COX6B2 were performed. The affinity constant is 11.7 nM.