

For Research Use Only

SLC22A4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83192-2-PBS

Basic Information

Catalog Number: 83192-2-PBS	GenBank Accession Number: BC028313	Purification Method: Protein A purification
Size: 100ug , Concentration: 1 mg/ml by Nanodrop;	GeneID (NCBI): 6583	CloneNo.: 240064F6
Source: Rabbit	UNIPROT ID: Q9H015	
Isotype: IgG	Full Name: solute carrier family 22 (organic cation/ergothioneine transporter), member 4	
Immunogen Catalog Number: AG3974	Calculated MW: 551 aa, 62 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse, human

Product Information

83192-2-PBS targets SLC22A4 as part of a matched antibody pair:

MP00292-2: 83192-2-PBS capture and 83192-4-PBS detection (validated in Cytometric bead array)

Unconjugated rabbit recombinant monoclonal antibody in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

This conjugation ready format makes antibodies ideal for use in many applications including: ELISAs, multiplex assays requiring matched pairs, mass cytometry, and multiplex imaging applications. Antibody use should be optimized by the end user for each application and assay.

Storage


Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00292-2, SLC22A4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83192-2-PBS. Detection antibody: 83192-4-PBS. Standard: Ag3974. Range: 1.25-40 ng/mL.