

For Research Use Only

# PROCA1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83193-3-PBS

## Basic Information

<b>Catalog Number:</b> 83193-3-PBS	<b>GenBank Accession Number:</b> BC029574	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 100ug , Concentration: 1 mg/ml by Nanodrop;	<b>GeneID (NCBI):</b> 147011	<b>CloneNo.:</b> 230538D11
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q8NCQ7	
<b>Isotype:</b> IgG	<b>Full Name:</b> proline-rich cyclin A1-interacting protein	
<b>Immunogen Catalog Number:</b> AG10851	<b>Calculated MW:</b> 336 aa, 38 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human


## Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve ofMP00294-2, PROCA1 Recombinant Matched Antibody Pair - BSA and Azide free. Capture antibody: 83193-3-PBS. Detection antibody: 83193-4-PBS. Standard: Ag10851. Range: 5-40 ng/mL.